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
*Array Tomography in Electron Microscopy is one way to acquire data for the 3D reconstruction of structures. Thereby, serial sections are produced and collected on a substrate for subsequent imaging. Depending on the size or volume of a targeting region more or less sections are required for reconstruction.*

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


*Array Tomography allows imaging of the same section at different pixel resolutions. Images taken of low resolution (700 nm/pixel), intermediate resolution (30 nm/pixel), and high resolution (3 nm/pixel).*

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*In comparison to other 3D-techniques, such as Serial Block-Face Scanning EM (SBFSEM) and Focused Ion Beam (FIB), Array Tomography is non-destructive. Collected sections can be stored, re-imaged at difference regions and resolutions, re-stained with heavy metals to enhance contrast or with fluorescent markers for correlative light and electron microscopy (CLEM).*

*Array Tomography is used in various research fields not only for reconstruction of biological structures but also for complex materials.*

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#### ► Non-Destructive

*Sections are stored on substrate and can be kept for years.*

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#### ► Sample Library

*Stored sections can be re-evaluated at different locations and resolutions.*

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


#### ► Save Scan Time and Data Storage



*No need for highest resolution or whole area scans due to re-imaging at various resolutions and regions.*

► Large Volume Reconstruction with ATUMtome



▪ (a) Image of cortical column neurones in the cerebral cortex. (b,c) About  $1000\mu\text{m}^3$  around an apical dendrite of one pyramidal cell. This region contains around 675 synapses, 530 axons, and 90 dendrites with highly non-random connectivity.

▪ Images of tissue sections reconstructed through 3D imaging and array tomography.


*Images supplied courtesy: Jeff Lichtman, Harvard University, Boston, USA*

► Small Volume Reconstruction with ASH

▪ Root tip of *Arabidopsis thaliana*. (a) LM and EM image of section ribbons on ITO substrate. (b) Reconstruction of root tip vacuoles. Large vacuole in yellow wire mesh, small vacuoles in light blue. Volume rendering performed with UCSF Chimera.

*Images supplied courtesy: Irene Wacker and Rasmus Schröder, BioQuant, Heidelberg, Germany.*

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*ATUMtome – A fully automated system for large volume reconstructions*

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*In classical ultramicrotomy, serial sections are cut and float on the water surface in the diamond knife trough and are picked up manually.*

*The ATUMtome is a unique system for automated collection of sections on a continuous tape. Thereby, the tape moves through the knife trough and automatically collects the sections from the water surface in sequential order. After section collection the tape is cut in strips and is mounted on substrates such as silicon wafers or glass slides. Data acquisition is done with a scanning electron microscope.*

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*The ATUMtome in conjunction with the PowerTome 3D can produce many thousands of sections uninterrupted thanks to the 1mm advance of the PT3D. This means no loss of sections which is crucial for array tomography.*

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### ► **Fully Automated**

*The tape collection and the cutting process are synchronized. Tape and cutting speed can be adjusted individually.*

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### ► **Unattended Section Collection**

*Fully automated section collection process. Once set and stable, section pick-up by tape without the need of human intervention.*

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### ► **Easy Handling**

*Intuitive and simple section collection mechanism. No long training period is required.*

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### ► **Keeps Microscope Free**

*Sectioning process does not interfere with image acquisition. In SBFSEM or FIB sectioning is executed in the microscope.*

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*Controlled Environment Considerations*

Stabile Temperatur: ± 1°C /h <i>Stable temperature: ± 1°C /h</i>	Abschirmung gegen Windströmungen, Lüftungsanlagen <i>Shield from wind currents, ventilation systems</i>
Rel. Luftfeuchtigkeit: ~50% <i>Rel. Humidity: ~50%</i>	Verhindern der Probenblock-Aufladung beim Schneiden (DiATOME Ionisator) <i>Prevent sample block charging during sectioning (DiATOME Ionizer)</i>

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*DiATOME trim 90, for orthogonal trimming for serial sectioning*

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*DiATOME ultra 35° 4.0mm, for sectioning of large sample blocks*

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### ASH2 – The Advanced Substrate Holder of the 2<sup>nd</sup> Generation

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*The ASH2 is a multi-axis manipulation tool designed to precisely position substrates in a knife trough and is capable of collecting tens to hundreds of serial sections for array tomography and correlative 3D workflows.*

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*A substrate is clamped to the ASH2 and placed in the water of the knife trough. After sectioning, the water is slowly removed and the section ribbons adhere to the substrate.*

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*The easy setup, alignment and handling of the ASH2 allows an efficient collection of sections on a desired substrate even for small section series. Subsequent data acquisition is executed at the light (LM) or scanning electron microscope (SEM).*

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*The main application of the ASH2 is the take-up of a few hundred sections for reconstruction of small sample volumes.*

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### ► **Fast Set Up**

*Quickly mounts onto the knife stage, converting a standard ultramicrotome to an array tomography system.*

### ► **No Mounting Rail**

*No need for a mounting rail and special table like for the precursor ASH-100.*

### ► **Fits Any Ultramicrotome**

*Integrates with current RMC and Leica ultramicrotomes including UC6 and UC7.*

### ► **No Limitation to Knives**

*Compatible with commercially available diamond knives for room temperature. Usage with a variety of standard diamond knives as well as knives with larger troughs for large substrates.*

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*A simple solution for manual collecting of section ribbons is the ultra ats knife. This DiATOME diamond knife is characterized by a particularly large water trough that can be mounted with a substrate. The ribbons are manually maneuvered over the substrate. By lowering the water level through a hole in the bottom, the sections adhere to the substrate and are ready for further processing.*

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**Tech@ScienceServices.de**

# Any further Questions?

*If you have any further questions regarding a product or the field of microtomy, please do not hesitate to contact us. We also offer on-site demonstrations.*

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**Science Services GmbH**